

# Electrospray and Atmospheric Pressure Chemical Ionization Quadrupole MS for Small Molecules

Jonathan A. Karty, Ph.D.

<http://msf.chem.indiana.edu>



# LC-MS Picture



# General Sample Guidelines

- **Purify analyte if possible**
  - Analyte should be 1-100 mg/L in concentration
    - Colored samples should be translucent
  - *Sample must have no particulates or precipitates*
- **Use only volatile solvents/buffers**
  - MeOH, H<sub>2</sub>O, acetone, CH<sub>3</sub>CN, THF, CHCl<sub>3</sub>, etc.
  - HOAc, HCOOH, NH<sub>4</sub>OAc, NH<sub>3</sub>, etc. (weak, volatile)
  - Ionic strength < 20 mM is best (e.g. 0.1% v/v HOAc)
  - *1 – 5 μM NaOAc for acid/base labile samples*
- **Need at least 50 μL for loop injection**
- **If you need non-polar solvent for APCI, see Dr. Karty first**

# Agilent 6130 System Description

- Inlet is an Agilent 1200 HPLC System
  - Binary gradient (14 solvent bottles, only 2 can be used in an analysis (A1 or A2 with B1 or B2, not A1 and A2)
    - 50-2,500  $\mu\text{L}/\text{min}$  flow rates
  - Autosampler 0.5-100  $\mu\text{L}$  injection volume
- 6-column selection valve (loop or 1 of 5 HPLC columns)
  - $\text{C}_{18}$  reversed phase is default column ( $\text{H}_2\text{O}-\text{CH}_3\text{OH}$  solvents)
- Agilent 1200 DAD detector
  - 190-800 nm UV-VIS detector
  - Can record entire UV-VIS spectrum throughout a run
- Electrospray Ionization-Atmospheric Pressure Chemical Ionization/Quadrupole (ESI/APCI-Q) mass spectrometer
  - 50 – 3,000  $m/z$  range
  - Can make ions by ESI, APCI, or both
  - Can alternate positive and negative ion modes during a run

# Air Sensitive MS

- Source on 6130 MS is flushed constantly with gas from a liquid nitrogen dewar
  - Enables walk-up air sensitive ESI/APCI MS
- A syringe pump with air-tight “Sample-Lock” syringe can be used for air-sensitive analyses
  - Requires special training from Dr. Karty
  - A separate syringe is needed for each sample
- A “dummy” sample must still be placed in the autosampler
- Remember to keep concentrations around 20 mg/L
  - Colored solutions should be translucent


# HPLC-MS vs. Loop Injection

- MSF uses the 1200 autosampler for nearly all samples
  - Must inject something, even if syringe pump is used
- Sample can go through an empty loop
  - Solution must contain only sample and volatile components
  - Allows 2.5 minute analysis time
- Sample can also be separated by HPLC
  - Up to 5 columns are available
  - Currently only C<sub>18</sub> reversed phase available for general use

# Electrospray Ionization (ESI)

- Dilute solution of analyte (1 mg/L) infused through a fine needle in a high electric field
- Very small, highly charged droplets are created
- Solvent evaporates, droplets split and/or ions ejected to lower charge/area ratio
  - Warm nebulizing gas accelerates drying
  - Basic sites are eventually protonated  $(M+H)^+$
  - Ketones, esters, and alcohols can sodiate  $(M+Na)^+$
  - Acidic sites can be deprotonated  $(M-H)^-$
  - Chlorinated solvents can add a chloride  $(M+Cl)^-$
- Free ions are directed into the vacuum chamber

# ESI Picture



[http://newobjective.com/images/electro/spraytip\\_bw.jpg](http://newobjective.com/images/electro/spraytip_bw.jpg)



# Characteristics of ESI Ions

- ESI is a thermal process (1 atm in source)
  - Little fragmentation due to ionization (cf EI)
- Solution-phase ions are often preserved
  - e.g. organometallic salts
- ESI ions are generated by ion transfer
  - $(M+H)^+$ ,  $(M+Na)^+$ ,  $(M-H)^-$ , or  $(M+Cl)^-$ , rarely  $M^{+\bullet}$  or  $M^{-\bullet}$
- ESI often generates multiply charged ions
  - $(M+2H)^{2+}$  or  $(M+10H)^{10+}$
  - Most ions are 500-1500 m/z
  - ESI spectrum x-axis must be mass/charge (m/z or Th, not amu or Da)

# Advantages of ESI

- Gentle ionization process
  - High chance of observing intact molecular ion
  - Very labile analytes can be ionized
  - Non-covalent complexes can be studied by ESI
- Molecule need not be volatile
  - Proteins/peptides easily analyzed by ESI
  - Salts can be analyzed by ESI
- Easily coupled with HPLC
- Both positive and negative ions can be generated by the same source

# ESI Disadvantages

- Analyte must have an acidic or basic site
  - Hydrocarbons and steroids not readily ionized by ESI
- Analyte must be soluble in polar, volatile solvent
- ESI is less efficient than other sources
  - Most ions don't make it into the vacuum system
- ESI is very sensitive to contaminants
  - Solvent clusters can dominate spectra
- Distribution of multiple charge states can make spectra of mixtures hard to interpret
  - e.g. polymer mass spectra


# What Samples Can Be Analyzed by ESI-QMS?

- Non-volatile organic molecules
  - Should be soluble in polar, volatile solvents
  - Molecule must be ionizable (MS detects ions)
    - R-NH<sub>2</sub>, R-CO<sub>2</sub>H, R-HSO<sub>3</sub>, R-OH, R-H<sub>2</sub>PO<sub>3</sub> work best for ESI
    - Aromatics, ketones, protected heteroatoms work well by APCI
- Organometallic complexes
  - Organometallic salts work especially well
  - Source is flooded with N<sub>2</sub> gas from a dewar
- HPLC-MS
  - With appropriate mobile phases

# What Samples Are Inappropriate for ESI-QMS analysis?

- Samples in non-polar or non-volatile solvent:
  - Hexane, benzene,  $\text{CH}_2\text{Cl}_2$ , DMSO, etc.
  - Run these in APCI or mixed mode
- Buffer systems incompatible with ESI
  - 6M urea, 10% glycerol, 0.1 M  $\text{NaH}_2\text{PO}_4$ , TBAF, (involatile)
  - Strong acid/base solutions (e.g. 10 mM HCl or NaOH too conductive)
  - detergents
- Molecules that have no ionizable groups

# ESI Example



Chemical Formula:  $C_{77}H_{136}NaO_{13}Si_4^+$

Exact Mass: 1403.8950

Chemical Formula:  $C_{77}H_{134}NaO_{12}Si_4^+$

Exact Mass: 1385.8845




# Atmospheric Pressure Chemical Ionization (APCI)

- APCI uses a corona discharge to generate acidic solvent cations from a vapor
  - Ionizing reagent often  $\text{CHCl}_3^+$  or  $\text{CH}_3\text{OH}^+$
- These solvent cations can protonate hydrophobic species not amenable to ESI
  - APCI can be done from hexane or THF
  - Often used to study lipids and steroids
  - In MSF, completely protected macrocycles are routinely studied by APCI
- APCI is a little harsher than ESI
  - Labile molecules might fragment



# APCI Diagram

ATMOSPHERIC PRESSURE




<http://imaisd.usc.es/riaidt/masas/imagenes/apci1.jpg>

# APCI Example


-2-43

j-2-43b 4 (0.204) AM (Cen,4, 80.00, Ar,10000.0,0.00,0.00); Cm (3:6)

2: TOF MS AP+  
2.39e6




j-2-43b 4




Chemical Formula:  $\text{C}_{22}\text{H}_{27}\text{F}_2\text{O}_3^+$   
Exact Mass: 377.1923

1: TOF MS ES+  
4.68e6



# Agilent 6130 Multi-mode Source




# Cluster and Background Ions

- The low  $m/z$  (<300) region of the mass spectrum is often dominated by solvent cluster ions
  - These are distinguished by low mass defects due to  $\text{Na}^+$  and  $\text{K}^+$  incorporation
  - Masses are  $\text{XXX}.0$  or  $\text{XXX}.9$
- Compounds from previous samples can coat the source leading to persistent “background” ions
  - e.g. 242.1 for tetrabutylammonium
  - Might see large ion from previous person’s run
- To confirm an ion is “background” try injecting a blank of ESI solution with your solvent system

# About the Fragmentor...

- There is a region in the source where ion-nitrogen collisions can cause fragmentation
  - This occurs between the end of the heated capillary and the first skimmer cone
- The voltage difference between capillary and skimmer is called the “Fragmentor” voltage
- Fragmentor voltage has 2 helpful main effects
  - It increases the number of ions that make it to the quad
  - It can knock non-covalent attached neutrals
- However, a high fragmentor value can cause labile groups to dissociate
  - Phospho-diesters can be rather vulnerable to this

# Source Diagram



# A Word About Quantification

- LC-MS is a quantitative technique
- Use only 1 m/z when quantifying a compound
- Intensity is proportional to concentration
  - $I \propto [X]$
  - $\alpha$  is unique to each compound
  - **The more two compounds differ chemically, the more careful one must be when comparing their intensities**
- Ideally a calibration curve is constructed using multiple solutions of pure analyte at varying concentrations
- A calibration lecture will be given next month if there is sufficient interest

# Hands-on Training

- Training starts soon, contact Jon or Angie
- Groups of no more than three
- One hour or so to complete
- No charge for first session
- After training, students must demonstrate competency by running their own samples prior to being granted after-hours access